


HOLO MODULE

Holographic device for measurement of semiconductor photoelectric parameters




HOLO-module is a novel device for contactless investigation and optical diagnostics of carrier generation, transport and relaxation processes in bulk semiconductor crystals, epilayers, and heterostructures. It explores a **four-wave mixing configuration** of the well-known “**pump-probe**” **technique**. The operation is based on the recording of a refractive index spatial modulation by interference field of the two picosecond laser pulses

and subsequent probing of the modulation decay by the delayed probe beam. **HOLO device** allows measurement of *carrier generation rate*, photoinduced *carrier concentration*, *bipolar diffusion coefficient*, *carrier recombination time* and *surface recombination rate*. The device also allows planar monitoring (mapping) of diffraction efficiency and extracting the homogeneity of defect distribution in a wafer.

BASICS OF OPERATION

Two coherent laser pulses overlap in the sample at angle 2Θ and create an interference pattern with period $\Lambda = \lambda/2\sin(\Theta/2)$. The interband transitions or impurity-assisted carrier generation creates a spatially modulated nonequilibrium carrier distribution $N = N_0(1 + \cos 2\pi x/\Lambda)$, which, in turn, is followed by a periodic refractive index modulation $\Delta n(x,t) \sim \Delta N(x,t)$. In this way, a dynamic free carrier grating is recorded in a semiconductor. The grating is probed by a delayed third beam at the wavelength far from the bandgap, and thus is partially diffracted on the grating. The energy of the diffracted I_1 and transmitted I_T parts of the probe are recorded in the far field of diffraction, and the instantaneous diffraction efficiency of



the grating $\eta = I_1/I_T$ is calculated. The diffraction efficiency, being nonlinearly dependent on the modulated part of the refractive index, $\eta \sim n^2$, sensitively reflects the changes of ΔN , thus providing carrier relaxation kinetics. Therefore, the HOLO-module allows optical studies of the carrier transport and determination of number of photoelectrical parameters, important for evaluation of semiconductor technology.

FEATURES

- Grating recording wavelength 266, 355, 532 or 1064 nm
- Probe beam wavelength 532 or 1064 nm
- Excitation energy flux density 0.1–10 mJ/cm²
- Probe beam delay range 1500 ps
- Laser pulse duration range 5–50 ps
- Dynamic grating period range 2–20 μm
- Minimal probe beam diameter 100 μm

RANGE OF MATERIAL PARAMETERS

- Diffusion coefficient D – 0.1–50 cm²/s
- Carrier lifetime τ_R – 0.1–10 ns
- Surface recombination velocity S – 10⁴–10⁶ cm/s

EXAMPLES OF PARAMETERS DETERMINATION

DETERMINATION OF DIFFUSION COEFFICIENT D AND CARRIER LIFETIME τ_R

Diffraction efficiency η as a function of probe delay time Δt is called the grating decay kinetics, $\eta(\Delta t)$. The kinetics is measured for a selected grating period Λ , and the grating decay time τ_G (the time interval in which the diffraction efficiency decreases by e^2 times) is determined. In this time interval, the carrier modulation decays by e times due to carrier recombination and diffusion according to relationship $1/\tau_G = 1/\tau_R + 4\pi^2 D/\Lambda^2$. Thus, measuring a set of diffraction efficiency kinetics at different periods Λ allows to plot a function $1/\tau_G = f(\Lambda^2)$ (so called an angular characteristic) and extract the D and τ_R values from a linear fit of the latter plot.




Figure a). The grating decay kinetics at various grating periods, measured in heavily doped GaAs:C double heterostructures, manufactured for HBT.




Figure b). The angular characteristics for two differently doped GaAs:C double heterostructures. Linear fit of the latter characteristics provides the D and minority carrier mobility values in heavily doped GaAs:C layers.

DETERMINATION OF SURFACE RECOMBINATION VELOCITY S

At band-gap excitation, the carriers are excited in a very thin surface layer ($\sim 0.1 - 0.2 \mu\text{m}$), therefore surface recombination may take place until car-

rier diffuse away from the sample surface. The decay kinetics of a surface grating allows determination of a surface recombination velocity by numerical fit-


ting of the non-exponential grating decay with a solution of two-dimensional continuity equation.

DETERMINATION OF CARRIER GENERATION RATE

Measurement of diffraction efficiency vs. excitation density (so called exposure characteristic) provides an important information about carrier generation and recombination. A power law dependence of diffraction efficiency vs. excitation is characterized by its slope $\gamma = \Delta \log(\eta) / \Delta \log(I_0)$. γ value may vary from $\gamma = 2$ for linear interband carrier generation to $\gamma = 4$ for two-phonon car-

rier generation. The nonlinear recombination may diminish γ value causing a saturation of the exposure characteristics.

Figure c). Exposure characteristics in 50 nm thick InGaN layers with different In content. The saturation of the characteristics at $\sim 1 \text{ mJ/cm}^2$ indicates the threshold of stimulated emission in InGaN layers.



LASER LIGHT SOURCE

Ekspla **PL2140 series** picosecond lasers are highly recommended.

- High pulse energy at 266, 355, 532, 1064 nm
- Excellent pulse energy (<1,5 %) and duration (<1,0 %) stability
- All solid-state mode-locking

Requests for custom made products are welcome.

Contact person for more detail information:

Kęstutis Jarašiūnas,
e-mail: kestutis.jarasiunas@ff.vu.lt



E K S P L A Ltd.
Lasers and Laser Systems Div.
Savanoriu av. 231
02300 Vilnius – 53
L I T H U A N I A

Ph.: +370 5 2649629
Fax: +370 5 2641809
sales@ekspla.com
www.ekspla.com

SALES AND SERVICING REPRESENTATIVES

USA
ALTOS Inc.
Ph.: +1-949-216 9012
Fax: +1-909-363 8637
sales@altos-inc.com

ISRAEL
BEN MOSHE
Ph.: +972-3-670 0007
Fax: +972-3-672 7319
bms@shani.net

JAPAN
TOKYO INSTRUMENTS
Ph.: +81-3-368 647 11
Fax: +81-3-368 608 31
sales@tokyoinst.co.jp

GREAT BRITAIN
PHOTONICS SOLUTIONS plc.
Ph.: +44-0131-664 8122
Fax: +44-0131-664 8144
sales@psplc.com

GERMANY
TOPAG LASERTECHNIK GmbH
Ph.: +49-6151-425 978
Fax: +49-6151-425 988
info@topag.de

KOREA
WOORYANG Co.
Ph.: +82-2-525 3816
Fax: +82-2-586 6721
wooyang@coslaser.com

FRANCE
BFI OPTILAS
Ph.: +33-(0)1-607 95936
Fax: +33-(0)1-607 98901
dadi.wang@bfioptilas.com

AUSTRALIA
FIBERBYTE
Ph.: +08-8410 4955
Fax: +08-8410 4544
rshaw@fiberbyte.com

CHINA
BEIJING OPLAN S&T Co., Ltd.
Tel: +86-10-62623871
Fax: +86-10-62616041
mxh@oplanchina.com

INDIA
ANATECH INSTRUMENTS
Ph.: +91-22-26730463/26730464
Fax: +91-22-26730465/26932893
anatech@bol.net.in

TAIWAN
SUPERBIN Ltd.
Ph.: +886-2-270 136 26
Fax: +886-2-270 135 31
davidyang@superbin.com.tw